



# Trusted

---

## Vapomatic<sup>TM</sup>

### Anesthetic Vaporizers

So dependable, they wear our family name.



**AM BICKFORD**  
SINCE 1974 • SAME NAME • SAME COMPANY inc.

800.795.3062 • [ambickford.com](http://ambickford.com)



800.795.3062

[ambickford.com](http://ambickford.com)  
[sales@ambickford.com](mailto:sales@ambickford.com)

**A.M. BICKFORD**  
SINCE 1974 • SAME NAME • SAME COMPANY inc.

### Features at a glance:




#### Manufactured in USA

Proudly made in the USA for more than 40 years.




#### Excellent Customer Service

Dedicated and dependable service, tailored to the needs of veterinarians.



#### On-time Shipping

Competitive rates.  
Orders ship daily.  
Local pickup also available.




#### Large Easy-Read Dial

Large dial with easy to see numbers has "off" lock and calibration indents.




#### Tamperproof Adjustments

Dial seal must be broken for access.



#### Standard Fittings

Standard 23mm male and female tapers are individually replaceable.




#### Light Weight

Dry weight approximately 7.5 lb.



#### Temperature Compensated

Range of ambient temperatures: 60°F to 90°F.



#### Liquid Capacity

Liquid capacity (dry) is approximately 200 ml.

# Vapomatic™ - Anesthetic Vaporizers

[Click here to see Vapomatic™ Vaporizers online.](#)



**Isoflurane  
Vapomatic™ Vaporizer  
#20010**



**Sevoflurane  
Vapomatic™ Vaporizer  
#20020**

### Anesthesia with assurance.

Designed with the veterinarian in mind, the Vapomatic veterinary vaporizer is a superior performing, agent specific calibrated anesthetic vaporizer. The temperature compensated Vapomatic provides outstanding performance at the low flow rates commonly used for veterinary anesthesia. Available in Isoflurane or Sevoflurane

For more than 50 years, A.M. Bickford's veterinary anesthesia equipment and supplies have built a reputation of quality, safety, durability, and trust.



Visit [ambickford.com](http://ambickford.com)  
for information on our line  
of veterinary equipment & supplies.

800.795.3062

[sales@ambickford.com](mailto:sales@ambickford.com)



# Trusted Calibration

## Vapomatic™ Anesthetic Vaporizers

So dependable, they wear our family name.

Visit  
[ambickford.com](http://ambickford.com)



STEP 04  
Vapomatic™  
Shipped Back  
to customer

STEP 01



Contact A.M. Bickford

Call: 800-795-3062

Email: [sales@ambickford.com](mailto:sales@ambickford.com)

Use the online form.

**TIP!**  
FREE loaner  
available.

STEP 02

Ship  
Vapomatic™

A.M. Bickford, Inc.  
12318 Big Tree Rd.  
Wales Center,  
NY 14169-0201



**TIP!**  
Drain  
anesthetic  
before  
shipping.

STEP 03  
Vapomatic™  
Calibrated



**TIP!**  
Vaporizer  
should be  
calibrated  
annually.