



Safety

Omnicon F/AIR Anesthetic Gas Filter

So dependable, they wear our family name.



**A
M** **BICKFORD**
SINCE 1974 • SAME NAME • SAME COMPANY inc.

800.795.3062 • ambickford.com



800.795.3062

ambickford.com
sales@ambickford.com

A M BICKFORD
SINCE 1974 • SAME NAME • SAME COMPANY inc.

Features:

- Disposable
- Ready to use
- Minimal storage
- Compact size
- No Maintenance

Specifications:

- Adsorbent Surface Area:
200,000 m3/grams
- Particle Size, Average, U.S. sieve:
6 + 16
- Average Density:
48 grams/cc
- Void Volume:
40%
- Pressure Drop (psi):
<0.1
- Average Life:
12-15 hours
- Electrical Properties:
Non-conductive

Adsorption Capacity:

- *Carbon Tetrachloride Activity,
weight increase:
60%
- **Halogenated Gas Retentivity:
50 grams

WARNING: DO NOT USE WITH FLAMMABLE GASES

References: MEDICAL WORLD NEWS, Nov. 2, 1973,
Volume 14 No. 40, p. 24. Aorn Journal, Nov. 1973,
Volume 18, No. 5 pp.872-4.

* The steady state percentage increase in weight
of adsorbent after air - which has been saturated
with carbon tetrachloride at zero degrees C -
is passed through the adsorbent at 25 degrees C.

** The amount of halothane retained prior to break-
through (breakthrough denoting the appearance
of adsorbate in the effluent stream).

Omnicon F/AIR Anesthetic Gas Filter



[Click here to see Omnicon F/AIR Gas Filter online.](#)

Omnicon F/AIR Anesthetic Gas Filter #80120

Omnicon F/AIR is a sensible answer to anesthesia gas problems in the operating room. The F/AIR anesthesia gas filter was specifically designed to remove waste anesthesia gases such as ISOFLURANE, HALOTHANE, ENFLURANE, etc. from the operating room environment.

Dimensions:

- Height.....6 inches
- Weight of Active Components200 grams
- Diameter 3 1/4 inches
- Connector, outside diameter 7/8 inches male


Visit ambickford.com
for information on our line
of veterinary equipment & supplies.

800.795.3062

sales@ambickford.com

Pollution Control Using F/AIR Canister & Accessories

80010 PC-2A Non-Rebreathing System



Ordering Information:

- #80120 F/AIR Canister Unit, or in cases of 8
- #81915B Adapter Sleeve, 22mm, F x 15mm. F
- #80020 F/AIR Holder
- #08106 Endo, Tube Connector, 6mm
- #18007 Scavenging Tube 40" with 3 adapters
- #80005 PC-1 Pop-Off Valve
- #80010 PC-2A NRB System
- #81915 Adapter, 15mm. F x 22mm. M
- #11140 Blue Plastic Scavenging Tubing with 22mm cuffs at 6" intervals, specify length
- #80122 Set of 3 adapters: 81915, 81915B & 08106
- #80000 F/AIR Kit consisting of one each, 80120, 80020, and 18007



Visit ambickford.com
for information on our line
of veterinary equipment & supplies.

800.795.3062

sales@ambickford.com